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***Measurement and instrumentation***

**Ref :  
293 057**

**Français – p 1**

**English – p 17**

**Version : 9001**

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**NOTES**


- ① On – off switch
- ② Power-on indicator lamp
- ③ Voltage selector switch:  $3 \times 6 \text{ V} - 3 \times 9 \text{ V}$  RMS between the live and neutral
- ④ Frequency selector switch: 0.1 Hz - 1 Hz - 50 Hz - 150 Hz and “0.1 Hz” with manual displacement
- ⑤ Output terminals

## 1 Measuring the phase to ground and phase to phase voltages

In a three-phase installation there are three lives and one neutral wire.

The voltage between the phase and the neutral is called phase to ground voltage; its effective value is noted as:  $V$ .

The voltage between the two phases is called the line voltage; its effective value is noted as:  $U$ .

## 1.1 Objective

To check the relation between phase to ground voltage  $V$  and line voltage  $U$ :

$$U = V \times \sqrt{3} .$$

## 1.2 Equipment

- An extra-low voltage three-phase generator;
- Three digital voltmeters
- Six connecting wires.

## 1.3 Measurements

- Measure the phase to ground and line voltages; complete the table of measurements.

Example of results:

Rated capacity	Phase to ground voltages			Line voltages		
	$V_{1N}$	$V_{2N}$	$V_{3N}$	$U_{12}$	$U_{23}$	$U_{31}$
$3 \times 6 \text{ V}$						
$3 \times 9 \text{ V}$	$V_{1N} = 8.9 \text{ V}$	$V_{2N} = 8.9 \text{ V}$	$V_{3N} = 8.9 \text{ V}$	$U_{12} = 15.5 \text{ V}$	$U_{23} = 15.3 \text{ V}$	$U_{31} = 15.4 \text{ V}$

- Compare for each of the two rated capacities the values of the phase to ground voltages  $V$  then the line voltages  $U$ .
- Check the relation  $U = V \times \sqrt{3}$

## 2 Viewing the voltages using ESAO4 and Générís

### 2.1 Objective


- To view in three-phase sinusoidal mode: The voltages between the phase and the neutral  
The voltages between the phases.
- To verify that a three-phase voltage system is a system with three sinusoidal voltages having the same amplitude, same frequency and offset one from the other by a third of a period
- To verify that the line voltages are leading by  $\frac{\pi}{6}$  from the phase to ground voltages
- To construct Fresnel's diagram for the phase to ground and line voltages (see 4).

## 2.2 Equipment

- An extra-low voltage three-phase generator (9V - 1 Hz)
- A computer with Générís installed
- ESAO4 console
- Three voltmeter adaptors
- Six connecting wires.

## Phase to ground voltages


## 2.3 Set-up

## 2.4 Experiment

- Perform the assembly.
- Start the computer; click on the Générís icon (the indicator lamp on the console comes on).
- Select the Test and calibrate icon
  - Select: Channel: **Voltmeter /**  
Rated capacity:  25 V  
 Instantaneous
  - Adjust the zero
  - Click on OK
  - Select: Channel: **Voltmeter //**  
Rated capacity:  25 V  
 Instantaneous
  - Adjust the zero
  - Click on OK
  - Select: Channel: **Voltmeter ///**  
Rated capacity:  25 V  
 Instantaneous

- Adjust the zero
- Click on OK
- Select the Parameters icon
  - Validate: Acquisition mode: f(t) ▼  
 Acquisition time: 2 s ▼  
 Number of points: 501 ▼  
 Synchronisation
  - Select: Channels on the ordinate: Voltmeter/  
 Voltmeter//  
 Voltmeter///
  - Click on Synchrono...
  - Select: Channel: Voltmeter/ ▼  
 Level: 0  
 Synchronisation  Increasing
  - Click on OK
  - Click on OK
- Select the Launch icon
  - Select: Magnitude: v1  
 Names of the curves to be acquired: v1(t)  
 Magnitude: v2  
 Names of the curves to be acquired: v2(t)  
 Magnitude: v3  
 Names of the curves to be acquired: v3(t)
  - Click on Acquisition

## 2.5 Processing the curves

- Click on the Representation of curves icon
  - Select: Scales:  Automatic.
- Click on the Colours & Styles icon to modify the colour and the lines of the curves.

"The curves obtained seem to be the graphical representation of three sinusoidal functions with the same amplitude, same angular frequency and having a phase shift amongst them":

$$v_1(t) = \hat{U} \times \sin(\omega \cdot t)$$

$$v_2(t) = \hat{U} \times \sin\left(\omega \cdot t + \frac{2\pi}{3}\right)$$

$$v_3(t) = \hat{U} \times \sin\left(\omega \cdot t + \frac{4\pi}{3}\right)$$

- Click on the Pointer icon
  - Use the cross hair to measure:
    - The period of each phase to ground voltage (point to the "origin", drag the mouse to the "end" point by holding the left button down during the entire operation and read the value of the period in the Caption bar)
    - The offset from one phase to ground voltage to the other
    - The maximum value of each of the phase to ground voltages.

Use Zoom for a more accurate reading.

- Select Modelling then Function
  - Select:  
Curve to be modelled: u1(t) ▼  
 $v_1(t) = a \cdot \sin(\omega \cdot t)$   
Interval:                   from: 0 (for example)  
  to: 2 (for example)  
Name: V1(V) (for example) ;  
Parameters            a 12.8 (for example)  
    w 6.28 (for example)
  - Click on Plot
  - Adjust the curve by selecting      a 12.8 (for example)  
    w 6.28 (for example)
  - Click on Plot then on Yes

We obtain for example: Parameters   a  12.7  
   w  6.25  
   Standard deviation S = 319mV

- Click on Save then Close.

Carry out the same search for:

$$v_2(t) = a \cdot \sin(\omega t + b)$$


and

$$v_3(t) = a \cdot \sin(\omega t + b)$$

We obtain for example:

Parameters a  12.8  
 w  6.29  
 b  2.02  
 Standard deviation S =  
 201mV

Parameters a  12.9  
 w  6.27  
 b  4.14  
 Standard deviation S =  
 209mV

Check that  $v_1(t) + v_2(t) + v_3(t) = 0$  ( $\vec{V}_1 + \vec{V}_2 + \vec{V}_3 = \vec{0}$ ).

- Select Processing then Formula
  - Select: Curve to be processed: v1(t) ▼
 

Ordinate:	Name: v
	S.I. Units: V
Formula	v = : v1+v2+v3
Resulting curve:	v(t)
  - Click on Plot then Save and then Close
- Select Processing then Regression
  - Select: Curve to be processed: v(t) ▼
 


Processing interval	Start: 0
	End: 2
Equation	V(t)
  - Click on Plot

We read for example:

$$23.9E-3 t + 182E-3$$


Correlation coefficient:

0.2%

## Line voltages

### 2.6 Set-up


### 2.7 Experiment

- Perform the assembly.
- Select the Test and calibrate icon
  - Select: Channel: Voltmeter /  
 Rated capacity: Ⓞ 25 V  
 Ⓞ Instantaneous

- Adjust the zero
- Click on OK
- Select: Channel: Voltmeter//  
Rated capacity:  25 V  
 Instantaneous
- Adjust the zero
- Click on OK
- Select: Channel: Voltmeter///  
Rated capacity:  25 V  
 Instantaneous
- Adjust the zero
- Click on OK
- Select the Parameters icon
  - Validate: Acquisition mode, f(t) ▼  
Acquisition time: 2 s ▼  
Number of points: 501 ▼  
 Synchronisation
  - Select: Channels on the ordinate: Voltmeter/  
Voltmeter//  
Voltmeter///
  - Click on Synchro...
  - Select: Channel: Voltmeter/ ▼  
Level: 0  
Synchronisation   
 Increasing
  - Click on OK
  - Click on OK
- Select the Launch icon
  - Select: Magnitude: u12  
Names of the curves to be acquired: u12(t)  
Magnitude: u23  
Names of the curves to be acquired: u23(t)  
Magnitude: u31  
Names of the curves to be acquired: u31(t)

- Click on Acquisition

## 2.8 Processing the curves


- Click on the Representation of curves icon
  - Select: Scales:  Automatic.
  - Click on the Colours & Styles icon to modify the colour and the lines of the curves.

Represent  $u_{12}(\pi)$  ,  $u_{23}(\pi)$  and  $u_{31}(\pi)$ .

- Select Processing then Mapping
  - Select: Curve to be mapped: u12(t) ▼
 

Abscissa	S.I. Units:	pi
	Expression:	2*t
Ordinate:	Name:	u12
	S.I. Units:	V
	Expression	u12
Mapped curve:		u12(phi)
  - Click on Plot then Save and then Close

Do the same with the other curves  $u_{23}(t)$  and  $u_{31}(t)$ .


- Click on the Pointer icon
  - Use the cross hair to measure:
    - The period of each of the line voltages
    - The offset from one line voltage to the other
    - The maximum value of each of the line voltages.

Use Zoom for a more accurate reading.

## Phase to ground and line voltages

### 2.9 Set-up


### 2.10 Experiment

- Perform the assembly.
- Start the computer; click on the Générés icon (the indicator lamp on the console comes on).
- Select the Test and calibrate icon
  - Select: Channel: Voltmeter /  
Rated capacity:  25 V  
 Instantaneous
  - Adjust the zero
  - Click on OK
  - Select: Channel: Voltmeter //  
Rated capacity:  25 V  
 Instantaneous
  - Adjust the zero
  - Click on OK
  - Select: Channel: Voltmeter ///  
Rated capacity:  25 V  
 Instantaneous
  - Adjust the zero
  - Click on OK
- Select the Parameters icon
  - Validate      Acquisition mode: f(t) ▼  
                         Acquisition time: 2 ms ▼  
                         Number of points: 501 ▼  
                          Synchronisation
  - Select: Channels on the ordinate: Voltmeter /  
   Voltmeter //  
   Voltmeter ///


- Click on Synchro...
- Select: Channel: Voltmeter/ ▼  
Level: 0  
Synchronisation      Increasing
- Click on OK
- Click on OK
- Select the Launch icon
  - Select: Magnitude: v1  
Names of the curves to be acquired: v1(t)  
Magnitude: v2  
Names of the curves to be acquired: v2(t)  
Magnitude: u12  
Names of the curves to be acquired: u12(t)
  - Click on Acquisition

## 2.11 Processing the curves

- Click on the Representation of curves icon
  - Select: Scales:  Automatic.
- Click on the Colours & Styles icon to modify the colour and the lines of the curves.
- Click on the Pointer icon
  - Use the cross hair to measure:
    - The period for each of the phase to ground and line voltages
    - The phase shift in line voltages  $u_{12}$  in relation to phase to ground voltage  $v_1$
    - The maximum value of each of the voltages.

Use Zoom for a more accurate reading.

- Select Processing then Formula
  - Select: Curve to be processed: v1(t) ▼  
 Ordinate: Name: U12  
 S.I. Units: V  
 Formula U = : v1-v2  
 Resulting curve: U12(t)

- Click on Plot then on Save

Compare  $u_{12}(t)$  and  $U_{12}(t)$ .

### 3 Viewing the voltages on an oscilloscope

#### 3.1 Objective (see 2)


#### 3.2 Equipment

- An extra-low voltage three-phase generator
- A dual-trace oscilloscope
- Four connecting wires.

#### 3.3 Experiment


Work with a 150 Hz frequency to obtain a stable signal.

- ❖ View  $v_1$  and  $v_2$

- Press the INVERT Y2 button to view  $-v_2$
- Press the ADD button to view  $v_1 - v_2$
- Plot the oscillogram

- ❖ View  $u_{12}$

Compare the two signals  $v_1 - v_2$  and  $u_{12}$ .

To view the three voltages, use the two-way switch for phases 2 and 3.

## 4 Fresnel representation of voltages

We associate to every sinusoidal signal  $u(t)$  a Fresnel vector  $\vec{U}$  that rotates at an angular velocity  $\omega$ .


Represent the Fresnel vectors  $\vec{V}_1$ ,  $\vec{V}_2$  and  $\vec{V}_3$  associated to the three phase to ground voltages  $v_{1N}$ ,  $v_{2N}$  and  $v_{3N}$ .

Represent the Fresnel vector  $\vec{U}_{12}$  associated to the line voltage  $u_{12}$  knowing that:  $u_{12} = v_{1N} + v_{N2}$

That is  $\vec{U}_{12} = \vec{V}_1 - \vec{V}_2$

Graphically check the relation between the phase to ground and line voltage.

Such a system is called balanced.

## 5 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
**Rue Jacques Monod**  
**BP 1900**  
**27 019 EVREUX CEDEX FRANCE**  
**+33 (0)2 32 29 40 50**

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries