


# Magnetic Bumper Set

## ME-9885



Included Equipment	Part Number
Magnetic Bumper (2)	003-05806
Additional Equipment Recommended	
Force Accessory Bracket	CI-6545
Force Sensors	CI-6537 / CI-6746 / PS-2104

## Introduction

The Magnetic Bumper Set includes a pair of strong magnetic bumpers. These bumpers are ideal for measuring forces during elastic collisions between carts. Simply remove the hooks on the Force Sensors and replace with the magnetic bumpers. The bumper has a metal post which assures that it is mounted in the correct orientation.

## Suggested Activities

### Elastic Collisions

1. Remove the hooks from both Force Sensors
2. Screw a Magnetic Bumper into each of the sensors
3. Connect the Force Sensors to PASCO carts or Force Accessory Bracket



### Limited Warranty

For a description of the product warranty, see the PASCO catalog.

### Contacting Technical Support

Address: PASCO scientific  
10101 Foothills Blvd.  
Roseville, CA 95747-7100

Phone: (916) 786-3800  
(800)772-8700

### Copyright Notice

The PASCO scientific 012-09158 Magnetic Bumper Set Instruction Manual is copyrighted with all rights reserved.